CGUARD™

Carotid Embolic Prevention System

Table of Contents

No.	Language	
1	English	2
2	Česky	4
3	Deutsch	7
4	Dansk	10
5	Nederlands	12
6	Suomi	15
7	Français	18
8	Magyar	21
9	Italiano	23
10	Latviski	26
11	Lietuviškai	29
12	Norsk	31
13	Polish	34
14	Português	37
15	Русский язык	40
16	Srpski	43
17	Slovensky	45
18	Slovenščina	48
19	Español	50
20	Swedish	53

Instruction for Use English

1 DEVICE DESCRIPTION

The InspireMD CGuard™ Carotid Embolic Prevention System (EPS) is designed to deliver a self-expanding stent to the carotid arteries using a rapid exchange (Rx) delivery system. The self-expanding stent is constructed of a nickel titanium alloy (Nitinol) and is covered by a permanent protective mesh (MicroNet™). It is compressed and loaded into the Rx delivery system, delivered to the intended lesion site and then expanded by retraction of a protective sheath. The stent and mesh remain as a permanent vessel scaffolding implant. Upon deployment, the stent imparts an outward radial force on the arterial lumen to establish patency. The stents are available in the below size matrix.

Vessel diameter (mm)	Length (mm)				
	Diameter (mm)	20	30	40	60
4.8 – 5.7	6.0	CRX0620	CRX0630	CRX0640	CRX0660
5.6 – 6.5	7.0	CRX0720	CRX0730	CRX0740	CRX0760
6.4 – 7.3	8.0	CRX0820	CRX0830	CRX0840	CRX0860
7.2 – 8.1	9.0	CRX0920	CRX0930	CRX0940	CRX0960
8.0 – 9.0	10.0	CRX1020	CRX1030	CRX1040	CRX1060

Table 1 - CGuard™ Carotid EPS Size Matrix


Figure 1 - CGuard™ Carotid stent

The stent delivery system is a Rapid Exchange (Rx) system. The delivery system has a 135 cm working length and is compatible with 0.014" guidewires (195 cm length), 6F (ID > 2.24 mm) vascular sheath, or 8F (ID > 2.24 mm) guiding catheter, proximal protective devices with a 9F (ID > 2.12 mm) and 0.014" Distal Protective Devices.

The delivery system is comprised of an outer shaft and inner assembly. The Rx system is a pull-push platform available in one size for CGuard™ stents. A pictorial representation of the delivery system is presented in Figure 2.

Recommended accessories for use are: PT™ moderate support by Boston Scientific or its comparable, with at least 190 cm long and 0.014" guiding wire; 8F (ID > 2.24) Guiding catheter e.g. Mach 1™ by Boston Scientific or Cordis Vista Brite tip MPA 1 or its comparable, or guiding sheath 6F (ID > 2.24) Chariot by Boston Scientific or Neuron Max Vascular Long Sheath by Penumbra, or its comparable; Distal protection device such as Filter EZ™ by Boston Scientific, Spider FX™ by EV3 or their comparable; Proximal protective device such as a 9F (ID 2.12mm) MoMa ultra by Medtronic or its comparable.


Figure 2 - Rx Delivery system

2 HOW SUPPLIED

Sterile: This device is sterilized with ETO. Non-pyrogenic. Intended for one (1) use only.
Contents: One (1) CGuard™ Self-Expanding Carotid Stent with Rx Delivery System. The system is placed inside a hoop and secured to the tray for support.
Storage: Store in a dry, dark, cool place.

3 INDICATIONS

The CGuard™ EPS is indicated for:

Improving carotid luminal diameter in patients at high risk for adverse events from carotid endarterectomy who require carotid revascularization and meet both criteria outlined below:

- Patients with neurological symptoms and >50% stenosis of the common or internal carotid artery by either ultrasound or angiogram or patients without neurological symptoms and > 80% stenosis of the common or internal carotid artery by either ultrasound or angiogram.
- Patients having a vessel with reference diameters between 4.8 mm and 9.0 mm at the target lesion.

4 CONTRAINDICATIONS

The CGuard™ EPS is contraindicated for use in:

- Patients in whom anti-coagulant and / or anti-platelet therapy is contraindicated.
- Patients with severe vascular tortuosity or anatomy that would preclude the safe introduction of a guide catheter, sheath, or stent system.
- Patients with known hypersensitivity to nickel-titanium.
- Patients with uncorrected bleeding disorders.
- Lesions in the ostium of the common carotid artery.
- Pre-existing jaw claudication.

5 WARNINGS

5.1 General Warning

- Only physicians who have received appropriate training and are familiar with the principles, clinical applications, complications, side effects and hazards commonly associated with

carotid stent placement should use this device.

- As with any type of vascular implant, infection secondary to contamination of the stent may lead to thrombosis, pseudo-aneurysm, or rupture.
- Stenting across a major bifurcation may hinder or prevent future diagnostic or therapeutic procedures.
- When multiple stents are required, stent materials should be of similar composition.
- Do not use the product after the "Use By" date specified on the package.
- Do not use the product if the temperature indicator on inner pouch is black.
- Do not use contrast material while performing flush to the CGuard™ EPS.
- Maintain the patient's Activated Clotting Time (ACT) at > 250 seconds throughout system usage to prevent thrombus formation on the device.
- Maintain continuous flush while removing and reinserting devices on the guidewire.
- Perform all exchanges slowly to prevent air embolism or trauma to the artery.
- Caution should be used if pre-dilating the lesion without embolic protection as this may increase the risk of an adverse outcome.
- Implanting a stent may lead to dissection of the vessel distal and / or proximal to the stent and may cause acute closure of the vessel, requiring additional intervention (carotid endarterectomy, further dilatation, or placement of additional stents).
- The stent may cause a thrombus, distal embolization or may migrate from the site of implant down the arterial lumen.
- Appropriate sizing of the stent to the vessel is required to reduce the possibility of stent migration.
- In the event of thrombosis of the expanded stent, thrombolysis and PTA should be attempted.
- In the event of complications such as infection, pseudo-aneurysm or fistulization, surgical removal of the stent may be required.
- Overstretching of the artery may result in rupture and life-threatening bleeding.
- If a filter (embolic protection system) is used, allow for and maintain adequate distance between CGuard™ EPS and the filter to avoid potential engagement or entanglement. If filter engagement and/or entanglement or filter detachment occurs, surgical conversion or additional catheter based intervention may be required.
- If a filter (embolic protection system) is used, choose a system with at least 195 cm.

5.2 Device Use Warning

- This device is intended for single-use only. DO NOT use the product after the 'Use By' date noted on the packaging.
- Package contains one self-expanding carotid stent system compressed in Rx delivery system. Store at room temperature.
- DO NOT re-use. DO NOT re-sterilize, as this can compromise device performance and may increase the risk of cross-contamination due to inappropriate reprocessing.
- When using a proximal protective device the use of power injector system is not recommended. If power injector is used, CGuard™ EPS delivery system may move distally from its position.
- Ensure optimal positioning of the stent prior to deployment. Once deployment is initiated the stent cannot be repositioned or recaptured. Stent retrieval methods and the use of additional wires, snares, and/or forceps may result in additional trauma to the vessel or the vascular access site. Complications may result in bleeding, hematoma, pseudoaneurysm, stroke or death.
- Frequently observe the CGuard™ Carotid EPS under fluoroscopy during stent deployment.
- Overstretching of the artery may result in rupture and life-threatening bleeding.
- Implanting a stent may lead to dissection of the vessel distal to and/or proximal to the stent and may cause acute vessel closure requiring additional intervention such as carotid endarterectomy, placement of additional stents or further dilatation.
- Appropriate stent sizing is required to reduce the possibility of stent migration. The stent may migrate from the site of implant, embolize or cause a thrombus distally from the implant site down the arterial lumen.
- Maintain the patient's activated clotting time (ACT) at > 250 seconds throughout the stent delivery and implant procedure to prevent thrombus formation on the Rx delivery system. In the event of thrombosis of the expanded stent, thrombolysis and PTA may be attempted. If IIb/IIIa inhibitors are used, maintain the ACT at > 200 seconds.
- Maintain a continuous flush while removing and inserting devices over the guidewire of the embolic protection device. Perform all device exchanges slowly in order to prevent air embolism or trauma to the artery.
- Allow for and maintain an adequate distance between the embolic protection system and the stent delivery system and/or the deployed stent in order to prevent possible entanglement of the two systems. If entanglement should occur and cannot be corrected, surgical intervention should be considered.

6 PRECAUTIONS

6.1 General- Precautions

CAUTION: Venous access should be available during carotid stenting to manage bradycardia and/or hypotension by either pharmaceutical intervention or placement of a temporary pacemaker, if needed.

- The use of embolic protection devices or distal protection device is recommended when using CGuard™ EPS.

6.2 Stent Handling – Precautions

- Special care must be taken not to handle or in any way disrupt the stent on the delivery system. This is most important during catheter removal from packaging, placement over the guidewire and advancement through the hemostasis valve and guiding catheter.
- Carefully inspect the system to verify that the device has not been damaged. Do not use damaged product.
- Take care to avoid unnecessary handling which may kink or damage the delivery system. Keep the delivery system as straight as possible and the delivery handle stationary during deployment. Do not use if device is kinked.
- Ensure the hemostatic valve is open to ensure freedom of movement of the outer sheath during deployment.
- Do not expose the delivery system to organic solvents as structural integrity and/ or device function may be impaired.
- Do not attempt to deploy the stent from its delivery system while the system is not located in target lesion. If deployed, the stent cannot be retrieved back into the delivery system and the stent may become damaged.
- The stent and the delivery system are designed to perform as an integrated system and to be used only as designed.

6.3 Stent Placement – Precautions

- Ensure that the stent delivery system is fully flushed with heparinized saline prior to use. Do not use the delivery system if the flush is not observed exiting at the distal end of the catheter and through the locking valve.

- Use with bleedback control hemostatic valves is not recommended.
- The system is not compatible with guidewires or embolic protection devices larger than 0.014" wire (0.36 mm).
- The system must be used with a guiding catheter or sheath to maintain adequate support of the guidewire or embolic protection device throughout the procedure.
- If resistance is met during introduction or delivery of the stent delivery system, or during retraction of the outer sheath, the system should be carefully withdrawn and another system used.
- When the catheter is in the body, it should be manipulated only under fluoroscopy. Radiographic equipment that provides high quality images is needed.
- Do not attempt to pull a partially deployed stent back through the guiding catheter or sheath; dislodgement of the stent from the delivery system may occur.
- Withdraw the stent system carefully if resistance is encountered during movement through the sheath or during initial retraction of the sheath.
- The delivery system is not designed for use with power injection. Use of power injection may adversely affect device performance.
- When more than one stent is required to cover the lesion or if there are multiple lesions, the distal lesion should be stented first followed by stenting of the proximal lesion. Stenting in this order obviates the need to cross the proximal stent for placement of the distal stent and reduces the chance of dislodging stents that have already been placed.
- If overlap of sequential stents is necessary, the amount of overlap should be kept to a minimum (about 5 mm). In no instance should more than 2 stents overlap.

6.4 Post-Implant – Precautions

- Care must be exercised when crossing a newly deployed stent with other interventional devices to avoid disrupting the stent geometry and placement of the stent.
- In the event of thrombosis of the expanded stent, thrombolysis and PTA should be attempted.
- When using a proximal protective device it is required to re-sheath the tip (re-introduce the tip back into the delivery sheath of the CGuard™ EPS) to ensure safe system retrieval into the proximal protective device working channel.

7 MRI Safety Information

MR Conditional 

Non-clinical testing and MRI simulations were performed to evaluate the entire family of the CGuard Carotid Stent. Non-clinical testing demonstrated that the entire family of the CGuard Carotid Stent is MR Conditional. A patient with an implant from this family can be scanned safely in an MR system under the following conditions:

- Static magnetic field of 1.5-Tesla and 3-Tesla, only
 - Maximum spatial gradient magnetic field of 4,000-gauss/cm (40-T/m)
 - Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2-W/kg for 15 minutes of scanning (i.e., per pulse sequence) in the Normal Operating Mode
- Under the scan conditions defined, the CGuard Carotid Stent is expected to produce a maximum temperature rise of 5.0C after 15-minutes of continuous scanning (i.e., per pulse sequence). In non-clinical testing, the image artifact caused by the CGuard Carotid Stent extends approximately 5-mm from this device when imaged with a gradient echo pulse sequence and a 3-Tesla MR system. The lumen of the CGuard Carotid Stent cannot be visualized on the gradient echo or T1-weighted spin echo pulse sequences.

8 POTENTIAL ADVERSE EVENTS

Based on the literature and on clinical and commercial experience with carotid stents and embolic protection systems, the following list includes possible adverse events associated with these devices:

- Abrupt closure
- Acute myocardial infarction
- Allergic reaction (contrast medium; drug; stent or filter material)
- Amaurosis fugax
- Aneurysm or pseudoaneurysm in vessel or at vascular access site
- Angina Coronary ischemia
- Arrhythmia (including premature beats, bradycardia, atrial and/or ventricular tachycardia, atrial and/or ventricular fibrillation [VF])
- Asystole or bradycardia requiring placement of a temporary pacemaker
- Arteriovenous fistula
- Bleeding complications from anticoagulant or antiplatelet medication requiring transfusion or surgical intervention
- Cerebral edema
- Cerebral hemorrhage
- Cerebral ischemia
- Congestive heart failure (CHF)
- Death
- Detachment and/or implantation of a component of the system
- Dissection of blood vessel
- Distal embolic protection device thrombosis occlusion
- Emboli, distal (air, tissue, plaque, thrombotic material, stent)
- Emergent or urgent surgery (Carotid Endarterectomy [CEA])
- Emergent surgery to remove stent or distal embolic protection device
- Fever
- Hematoma at vascular access site, with or without surgical repair
- Hemorrhagic event, with or without transfusion
- Hyperperfusion syndrome
- Hypotension/Hypertension
- Infection, local or systemic including bacteremia or septicemia
- Ischemia/ infarction of tissue organ
- Pain (head/neck)/ severe unilateral headache
- Pain at catheter insertion site
- Renal failure/insufficiency secondary to contrast medium
- Restenosis of vessel in stented segment
- Seizure
- Stent distal embolic protection device entanglement/ damage
- Stent distal embolic protection device collapse or fracture
- Stent malapposition/migration
- Stent thrombosis occlusion
- Stroke / cerebrovascular accident (CVA) / transient ischemic attack (TIA)
- Total occlusion of the carotid artery
- Vascular thrombosis/occlusion at puncture site, treatment site, or remote site
- Vessel dissection, perforation or rupture
- Vessel spasm or recoil

9 STENT SIZE DETERMINATION

- Careful stent sizing is important to successful stenting. A minimum "interference" fit of 0.5 mm between the vessel and the stent is recommended in order to achieve optimum sizing and stent expansion of the self-expanding stent. For example, select a 6.0 mm stent to treat a 4.8 – 5.7 mm diameter vessel. Select a 7.0 mm stent to treat a 5.6 – 6.5 mm diameter vessel. The mean percentage of foreshortening for all stent sizes is less than 6%. The shortest stent length consistent with total lesion coverage is optimal. Should adequate coverage by one stent be impossible, a second stent may be used.
- The Rx delivery system has a 135 cm working length and is compatible with 0.014" guidewires. Sheaths and guiding catheters with a minimum inner diameter of 2.24 mm are recommended.

WARNING: The CGuard™ EPS is contraindicated for use with lesions in the ostium of the common carotid artery.

WARNING: Appropriate sizing of the stent to the vessel is required to reduce the possibility of stent migration.

10 MATERIAL REQUIRED

- 6F (ID> 2.24) vascular sheath or 8F (ID> 2.24) guiding catheter or sheath length should not interfere with stent Rx delivery system requirements.
- Optional balloon dilatation catheter.
- Any carotid EPD with a 195 cm and 0.014" guidewire.
- Two to three syringes (10 -20 cc).
- One 3 cc syringe.
- 500 cc heparinized saline solution (sterile).

CAUTION: CGuard™ EPS is not compatible with guidewires or embolic protection devices larger than 0.014" (0.36 mm).

11 PRE-PROCEDURE

Patient preparation and sterile precautions should be the same as for any angioplasty procedure. The placement of the carotid stent in a stenotic or obstructed carotid artery must be done in a procedure room with angiography capabilities. Angiography should be performed to map out the extent of the lesion and the collateral flow. Access vessels must be sufficiently patent or sufficiently recanalized to proceed with further intervention.

12 INSPECTION PRIOR TO USE

- Remove the CGuard™ EPS Self-Expanding Stent with Rx Delivery System from its protective packaging. Lay the device flat. Take care not to kink the shaft of the Rx delivery catheter system.

CAUTION: Carefully inspect the CGuard™ EPS to verify that the device has not been damaged in shipment. Do not use damaged equipment.

CAUTION: The Rx delivery system has an internal shaft. Take care to avoid unnecessary handling which may kink or damage the delivery system. Keep the Rx delivery system as straight as possible and the delivery handle stationary during deployment. Do not use if device is kinked.

- Inspect the Rx delivery system sheath to verify that it has not been damaged during Shipment and that the stent does not overlap the proximal marker.
- Ensure that the stent is fully covered by the sheath.

CAUTION: Special care must be taken not to handle or in any way disrupt the stent on the delivery system. This is especially important during delivery system removal from packaging, placement over the distal embolic protection device wire and advancement through a hemostatic valve and guiding catheter hub.

CAUTION: The stent on the Rx delivery system is intended to perform as a system. Do not remove the stent from the delivery system as removal may damage the stent. If removed, the stent cannot be put back on the Rx delivery system.

13 DELIVERY SYSTEM PREPARATION

CAUTION: Do not expose the CGuard™ EPS to organic solvents as structural integrity and/or function may be impaired.

- For device flushing use a 3 cc syringe filled with heparinized saline solution, maintain positive pressure until saline fluid drops are observed exiting the CGuard™ EPS at the distal end. This process may take 30 seconds.
- Alternatively: Fill an inflator with 20 ml heparinized saline. Connect the inflator to the CGuard™ EPS and maintain 5 atm pressure for 30 seconds, or until the system is flushed and fluid is observed exiting the CGuard™ EPS at the distal end.

CAUTION: Ensure that CGuard™ EPS is fully flushed with heparinized saline prior to use. Do not use the CGuard™ EPS if flush is not visible exiting at the distal end of the sheath.

CAUTION: Do not use contrast material while flushing.

Keep the device straight and flat to avoid kinking the shaft.

14 LESION PREPARATION

WARNING: Administer heparin dose sufficient to maintain an ACT of > 250 secs to prevent thrombus formation on the devices.

CAUTION: Venous access should be available during carotid stenting to manage bradycardia and/or hypotension by a pacemaker placement or pharmaceutical intervention, if needed.

CAUTION: The CGuard™ EPS must be used with a guiding catheter or sheath to maintain adequate support of the 0.014" guidewire or embolic protection device throughout the procedure.

CAUTION: Use of bleedback control hemostatic valves is not recommended.

CAUTION: When catheters are in the body they should only be manipulated under fluoroscopy. Radiographic equipment that provides high-quality images is needed.

WARNING: Maintain continuous flushing while removing and reinserting devices on the guidewire or embolic protection device. Perform all exchanges slowly to prevent air embolism or trauma to the artery.

- Place an embolic protection device or distal protective device
- If needed, pre-dilate the lesion with an appropriate size balloon dilatation catheter to a minimum of 3.0 mm after the distal protection device is in place beyond the lesion. *Note:* If no predilatation balloon is utilized, there must be a minimum luminal opening of 3.0 mm to enable retrieval of the tip of the CGuard™ delivery system.
- Maintain the embolic protection device stationary while withdrawing the balloon catheter.

15 PROCEDURE

- If lesion pre-dilatation has been performed, remove the balloon catheter and load the delivery system onto the 0.014" (0.36 mm) guide wire.

CAUTION: Using a proximal protection device (with minimum ID of 2.12 mm) along with a power injector for lesion visualization is not recommended, CGuard™ delivery system may move distally from its position during injection.

CAUTION: The Rx delivery system is not designed for use with a power injector. Use of a power injector may adversely affect device performance.

- Keep the device flat to avoid kinking the shaft.
 - Insert the Rx delivery system through the hemostatic valve adapter.
- CAUTION:** If resistance is met during Rx delivery system introduction, the system should be withdrawn and another system used.
- Advance the stent and Rx delivery system forward under fluoroscopic guidance to the lesion site.
- CAUTION:** Avoid tension in the Rx delivery system prior to deployment.

16 STENT DEPLOYMENT

WARNING: Ensure optimal positioning of the stent prior to deployment. Once deployment is initiated, the stent cannot be repositioned or recaptured. Stent retrieval methods (use of additional wires, snares and/or forceps) may result in additional trauma to the carotid vasculature and/or vascular access site. Complications may include death, stroke, bleeding, hematoma, or pseudoaneurysm.

CAUTION: Once stent placement has been initiated, do not attempt to pull a partially expanded stent back through the guiding catheter or sheath as dislodgement of the stent from the Rx delivery system may occur.

- Confirm the stent position angiographically prior to deployment. Adjust position if necessary.
- While holding the handle stationary with one hand, unlock the locking valve.

Note: Ensure that the Rx delivery system is straight and not coiled. Keep the inner shaft of the Rx delivery system stationary during deployment. Do not hold the outer sheath of the Rx delivery catheter during deployment. It must be free to move.

- Deploy the stent by slowly pulling back outer sheath. Continue to pull back until the stent is fully deployed.

Note: If significant resistance is encountered during pullback of the outer sheath and before stent release is initiated, re-lock valve and remove the system. Once deployment is initiated, the stent cannot be recovered by the sheath.

CAUTION: In the event of partial delivery of the stent as the result of the inability to fully deploy the stent, remove the entire Rx delivery system from the patient. This may result in damage to the vessel wall and may require surgical intervention.

- Rotate the valve to the locked position, carefully withdraw the distal tip of the Rx delivery system through the stent. Remove the Rx delivery system from the patient.
- Under fluoroscopy, confirm that the stent has been deployed at the target lesion.
- If additional stent-to-wall apposition is desired or to facilitate the use of other interventional devices, the stent can be post-dilated with a balloon dilatation catheter. Do not expand the stent beyond its unconstrained maximum diameter as stated on the label and in Table 1. Post-dilate as needed in accordance with the compliance chart accompanying the selected balloon catheter.

CAUTION: When more than one stent is required to cover the lesion or if there are multiple lesions, the distal lesion should be stented first followed by stenting of the proximal lesion. Stenting in this order avoid the need to cross the proximal stent in order to place the distal stent and reduces the chance of dislodging stents that have already been placed.

CAUTION: If overlap of sequential stents is necessary, the amount of overlap should be kept to a minimum (approximately 5 mm). In no instance should more than two (2) stents overlap.

CAUTION: Care must be exercised when crossing a newly stented lesion. Alternatively: stent with other interventional devices to avoid disrupting the stent geometry and placement of the stent.

WARNING: Overstretching of the artery may result in rupture and life-threatening bleeding.

17 POST STENT PLACEMENT

- Following stent deployment, carefully withdraw the delivery system until the tip crosses the lesion and then re-sheath (re-introduce the tip back into the delivery sheath of the CGuard™) the delivery system to ensure safe passage of the tip into the guiding catheter or the proximal protection device working channel. Then carefully withdraw the delivery system out of the patient body.

CAUTION: In the event the delivery system is not re-sheathed and closed, the tip may be caught in the working channel inlet and may be dislodged.

- Following stent placement, an angiogram should be performed to confirm vessel patency and percent stenosis remaining in the vessel lumen.

WARNING: The stent may cause a thrombus, distal embolization or may migrate from the site of implant down the arterial lumen. Appropriate sizing of the stent to the vessel is required to reduce the possibility of stent migration. In the event of thrombosis of the expanded stent, thrombolysis and PTA should be attempted.

- Upon completion of the angiogram, the embolic protection device should be removed in accordance with the instructions for use with that device.
- Patients should be put on an appropriate regimen of anticoagulants / antiplatelets.

WARNING: In the event of complications such as infection, pseudoaneurysm, or fistulization, surgical removal of the stent may be required.

18 WARRANTY/LIABILITY

The product and each component of its system have been designed, manufactured, tested and packaged with all reasonable care. The warnings contained in InspireMD instructions for use are expressly considered as an integral part of this provision. InspireMD warrants the product until the expiration date indicated on the same. The warranty is valid provided that the use of the product was consistent with the instructions for use. InspireMD disclaims any warranty of merchantability or fitness for a particular purpose of the product. InspireMD is not liable for any direct, indirect, incidental or consequential damages caused by the product. Except in the case of fraud or grave fault on InspireMD part, compensation of any damage to the buyer will not, in any event, be greater than the invoice price of the disputed products. The guarantee contained in this provision incorporates and substitutes the legal guarantees for defects and compliance, and excludes any other possible liability of InspireMD, however originating, from its product supplied.

These limitations of liability and warranty are not intended to contravene any mandatory provisions of law applicable. If any clause of the disclaimer is considered by a competent court to be invalid or to be in conflict with the applicable law, the remaining part of it shall not be affected and remain in full force and effect. The invalid clause shall be substituted by a valid clause which best reflects InspireMD legitimate interest in limiting its liability or warranty. No person has any authority to bind InspireMD to any warranty or liability regarding the product.

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CGUARD™
Carotid Embolic Prevention System

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InspireMD GmbH